


DR. AMIR CMD-CENTRUM CARE INSTRUCTIONS FOR THE DENTAL SPLINT



You just got a high-quality therapeutic dental splint from us. So it remains in good condition, please note the following care instructions:




THE RIGHT TEMPERATURE

Only clean the dental splint with **cold water**, as warm or hot water can cause deformation.

THE RIGHT BRUSH

For cleaning, use a separate, **soft toothbrush**.




THE RIGHT CLEANING PRODUCT

Use **dish soap or liquid soap** instead of toothpaste, as the latter will settle in the dental splint and will make it rough.

THE RIGHT STORAGE

After cleaning, store the damp or wet dental splint only in the **open box**. Do not close it immediately afterwards.



DO YOU HAVE ANY QUESTIONS ABOUT
YOUR DENTAL SPLINT?
FEEL FREE TO CONTACT US AT ANY TIME
BY PHONE OR VIA E-MAIL!

+49 511 / 544 24 44

verwaltung@cmd-hannover.de